

## anti- TP73 antibody

### Product Information

Catalog No.:	FNab10864
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Tumor protein p73
Synonyms:	p53 like transcription factor, p53 related protein, P73, TP73, tumor protein p73
Observed MW:	48-80kd
Uniprot ID :	O15350

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human esophagus using FNab10864(TP73 antibody) at dilution of 1:50



Hela cell lysates were subjected to SDS PAGE followed by western blot with FNab10864(TP73 antibody) at dilution of 1:1000