

## anti- PBX2 antibody

### Product Information

Catalog No.:	FNab10869
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	pre-B-cell leukemia homeobox 2
Synonyms:	G17, Homeobox protein PBX2, HOX12, PBX2, PBX2MHC, pre B cell leukemia homeobox 2, Protein G17
Observed MW:	45,60 kDa
Uniprot ID :	P40425

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Various cell lysates were subjected to SDS PAGE followed by western blot with FNab10869(PBX2 antibody) at dilution of 1:1000