

# anti- FITC antibody

#### **Product Information**

Catalog No.: FNab10886

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Goat

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.06% PC300 and 50% glycerol pH 7.3, -20°C for 12 months

(Avoid repeated freeze / thaw cycles.)

### **Background**

FITC (fluorescein isothiocyanate) is a fluorochrome dye that absorbs ultraviolet or blue light causing molecules to become excited and emit a visible yellow-green light. This emission ceases upon removal of the light causing the excitation. Fluorochrome labeling provides rapid, accurate localization of antigen-antibody interaction when one of the reactants is part of a cell, tissue or other biological structure. FITC is a commonly used marker for antibodies in immunofluorescent techniques since the conjugation of FITC to proteins is relatively easy and does not, in general, destroy the biological activity of the labeled protein. FITC is widely used as a hapten to label different proteins. Anti-FITC antibodies can be used for isolating FITC labeled cells by magnetic microbeads based sorting, as well as removal of background signal. The antibody reacts with FITC but not the carrier protein.

## Immunogen information

Immunogen: FITC conjugated to KLH
Synonyms: Fluorescein Isothiocyanate

**Application** 

Reactivity: FITC
Tested Application: ELISA

Recommended dilution: 0.5-1.0 µg/mL for elisa

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1