
anti- AMH antibody

Product Information

Catalog No.:	FNab10893
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9G2
IsoType:	IgG1
Storage:	PBS with 0.1% proclin300 pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	anti-Mullerian hormone
Synonyms:	AMH, anti Mullerian hormone, MIF, MIS
Uniprot ID :	P03971

Application

Reactivity:	Human
Tested Application:	Recommended for sandwich immunoassays in ELISA, CLIA, GICA.
Recommended dilution:	Capture antibody. Each laboratory should determine an optimum working titer for use in its particular application.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com