

anti- Phospho-NFKB1(S337) antibody

Product Information

Catalog No.:	FNab10905
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A phospho specific peptide corresponding to residues surrounding S337 of human NFKB1
Synonyms:	NFKB1, CVID12, EBP-1, KBF1, NF-kB1, NF-kappa-B, NF-kappaB, NFKB-p105, NFKB-p50, NFkappaB, p105, p50
Observed MW:	120 kd
Uniprot ID :	P19838

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Hela cells were subjected to SDS PAGE followed by western blot with FNab10905 at dilution of 1:1000



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab10905 at dilution of 1:1000