

anti- HDAC7 antibody

Product Information

Catalog No.: FNab10906

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: histone deacetylase 7

Synonyms: HD7, HD7A, HDAC7, HDAC7A, histone deacetylase 7, Histone

deacetylase 7A

Observed MW: 120 kd Uniprot ID: Q8WUI4

Application

Reactivity: Human, Mouse, Rat

Tested Application: WB

Recommended dilution: WB: 1:500-1:2000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Various lysates were subjected to SDS PAGE followed by western blot with FNab10906(HDAC7 Antibody) at dilution of 1:1000