

anti- PF4 antibody

Product Information

Catalog No.:	FNab10918
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Platelet factor 4
Synonyms:	C X C motif chemokine 4, CXCL4, Iroplact, Oncostatin A, PF 4, PF4, platelet factor 4, SCYB4
Observed MW:	10 kDa
Uniprot ID :	P02776

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Various lysates were subjected to SDS PAGE followed by western blot with FNab10918(PF4 Antibody) at dilution of 1:1000