
anti- Phospho-RB(S807/811) antibody

Product Information

Catalog No.:	FNab10920
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	A synthetic phosphorylated peptide around S807 & S811 of human RB
Synonyms:	OSRC, PPP1R130, RB, p105-Rb, pRb, pp110, RB1, Rb
Observed MW:	120 kDa
Uniprot ID :	P06400

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



HT-29 cells were subjected to SDS PAGE followed by western blot with FNab10920 at dilution of 1:1000. NIH/3T3 cells were treated by Nocodazole (100ng/mL) for 16 hours.