

---

## anti- IL27RA antibody

### Product Information

Catalog No.:	FNab10921
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Interleukin-27 receptor subunit alpha
Synonyms:	IL27RA CRL1, TCCR, WSX1, UNQ296/PRO336
Observed MW:	70 kDa
Uniprot ID :	Q6UWB1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Various lysates were subjected to SDS PAGE followed by western blot with FNab10921(IL27RA Antibody) at dilution of 1:1000