

anti- GTT1 antibody

Product Information

Catalog No.:	FNab10923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Glutathione hydrolase 1 proenzyme
Synonyms:	GGT1, CD224, D22S672, D22S732, GGT, GGT 1, GTG
Observed MW:	24kDa, 40-70kDa
Uniprot ID :	P19440

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Various lysates were subjected to SDS PAGE followed by western blot with FNab10923(GTT1 Antibody) at dilution of 1:1000