

anti- HIF2A antibody

Product Information

Catalog No.:	FNab10928
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Hypoxia-inducible factor 2-alpha
Synonyms:	EPAS1, ECTY4, HIF2A, HLF, MOP2, PASD2, bHLHe73, Endothelial PAS domain-containing protein 1
Observed MW:	120 kDa
Uniprot ID :	Q99814

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse brain tissue using FNab10928(HIF2A Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab10928 at dilution of 1:1000.