

## **anti- Phospho-eNOS(S1177) antibody**

### **Product Information**

Catalog No.:	FNab10931
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Immunogen:	A synthetic phosphorylated peptide around S1177 of human eNOS
Synonyms:	NOS3, ECNOS, eNOS
Observed MW:	133 kDa
Uniprot ID :	P29474

### **Application**

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	


### **Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



SKO3 cells were subjected to SDS PAGE followed by western blot with FNab10931 at dilution of 1:1000.