

## anti- CHRNE antibody

### Product Information

Catalog No.:	FNab10937
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Acetylcholine receptor subunit epsilon
Synonyms:	CHRNE, ACHRE, CMS1D, CMS1E, CMS2A, CMS4A, CMS4B, CMS4C, FCCMS, SCCMS
Observed MW:	45 kDa
Uniprot ID :	Q04844

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)

Various lysates were subjected to SDS PAGE followed by western blot with FNab10937( CHRNE Antibody) at dilution of 1:1000.

