

anti- SLC4A10 antibody

Product Information

Catalog No.:	FNab10942
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Solute carrier family 4 member 10
Synonyms:	SLC4A10, NBCn2, NCBE
Observed MW:	180 kDa
Uniprot ID :	Q6U841

Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution: WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



293T cells were subjected to SDS PAGE followed by western blot with FNab10942(SLC4A10 antibody) at dilution of 1:1000