

---

## anti- Phospho-STAT3 ( Y705 ) antibody

### Product Information

Catalog No.:	FNab10943
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	A synthetic phosphorylated peptide around Y705 of human Phospho-STAT3-Y705
Synonyms:	ADMIO, ADMIO1, APRF, HIES, STAT3, Stat3
Observed MW:	83 kDa,88 kDa
Uniprot ID :	P40763

### Application

Reactivity:	Human, Mouse
Tested Application:	WB
Recommended dilution:WB:	1:500-1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : ( 0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10943 at dilution of 1:500. HeLa cells were treated by IFN- $\alpha$  (100 ng/ml) at 37°C for 30 minutes after serum-starvation overnight.