

anti- PDGFRA antibody

Product Information

Catalog No.:	FNab10965
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Platelet-derived growth factor receptor alpha
Synonyms:	CD140A, CD140a antigen, PDGF R alpha, PDGF receptor alpha, PDGFR2, PDGFRA, PDGFR α , PDGFR- α , Rhe PDGFRA
Observed MW:	150 kDa
Uniprot ID :	P16234

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB :	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

