

anti- FLT3LG antibody

Product Information

Catalog No.:	FNab12000
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Fms-related tyrosine kinase 3 ligand
Synonyms:	Flt 3 ligand, Flt 3L, Flt3 L, FLT3 LG, Flt3 ligand, Flt3L, FLT3L_HUMAN, FLT3LG, Fms related tyrosine kinase 3 ligand, Fms-related tyrosine kinase 3 ligand, SL cytokine
Observed MW:	30kDa
Uniprot ID :	P49771

Application

Reactivity:	Human, Mouse
Tested Application:	WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Mouse intestine Lysates were subjected to SDS PAGE followed by western blot with FNab12000 at dilution of 1:600.