


Recombinant Mouse Prl3d1

Catalogue No.:	P2057
Species:	Mouse
Uniprot ID:	P18121
Expression Region:	30-224
Host:	E.Coli
Tags:	N-terminal His-IF2DI Tag
Molecular Weight:	43.0 kDa under reducing conditions
Purity:	Greater than 90% as determined by SDS-PAGE
Formulation:	Lyophilized from a 0.2 μ m filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH7.4
Reconstitution:	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to a concentration of 0.1-1mg/ml by gently pipetting 2-3 times, don't vortex.
Storage:	The lyophilized protein is stable at -20°C for up to 1 year, the protein solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle.

Synonyms:

DKFZp586G1822 antibody, DKFZp586P2220 antibody, FLJ41820 antibody, KIAA0618 antibody, MGC3792 antibody, Nuclear envelope pore membrane protein POM 121 antibody, Nuclear envelope pore membrane protein POM 121A antibody, Nuclear pore membrane protein 121 kDa antibody, Nucleoporin Nup121 antibody

SDS-PAGE:**Safety Note:**

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.