

Recombinant Human eIF2A

Catalogue No.:	P4511
Species:	Human
Uniprot ID:	P05198
Expression Region:	1-315
Host:	E.Coli
Tags:	N-terminal His Tag
Molecular Weight:	37.9 kDa under reducing conditions
Purity:	Greater than 95% as determined by SDS-PAGE
Formulation:	Lyophilized from a 0.2 μm filtered solution in 10 mM Hepes, 500 mM NaCl
	with 5% trehalose, pH7.4
Reconstitution:	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or
	buffer of interest to a concentration of 0.1-1mg/ml by gently pipetting 2-3
	times, don't vortex.
Storage:	The lyophilized protein is stable at -20°C for up to 1 year, the protein
	solution can be stored at 2-8°C for up to 1 week. For extended storage, it
	is recommended to further dilute in a buffer containing 0.1% BSA (make
	sure BSA will NOT interfere with your experiment) and store in working
	aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle.



Synonyms:65 kDa eukaryotic translation initiation factor 2A, CDA 02, CDA02, EIF 2,EIF 2A, eIF-2A, EIF2, eif2a, EIF2A_HUMAN, Eukaryotic translationinitiation factor 2A, Eukaryotic translation initiation factor 2A 65kDa,MST089, MSTP004, MSTP089

SDS-PAGE:

KDa 70 — 55 — 40 — 35 — 25 — 15 — 10 —

Safety Note: This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.