


Recombinant Human Lin28B

Catalogue No.:	P6070
Species:	Human
Uniprot ID:	Q6ZN17
Expression Region:	1-250
Host:	E.Coli
Tags:	N-terminal His Tag
Molecular Weight:	27.4 kDa
Purity:	Greater than 95% by SDS-PAGE gel analyses
Formulation:	Lyophilized from a 0.2 μ m filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH7.4
Reconstitution:	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 μ g/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Storage:	The lyophilized protein is stable at -20°C for up to 1 year, the protein solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle.

Synonyms:

CSDD 2, CSDD2, FLJ16517, Lin 28 homolog B (*C. elegans*), Lin 28.2,

Lin-28B, LIN28B, LN28B_HUMAN, Protein lin-28 homolog B

SDS-PAGE:**Safety Note:**

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.