


Recombinant Human ASB13

Catalogue No.:	P6305
Species:	Human
Uniprot ID:	Q8W XK3
Expression Region:	1-278
Host:	E.Coli
Tags:	N-terminal His Tag
Molecular Weight:	31.8 kDa under reducing conditions
Purity:	Greater than 95% as determined by SDS-PAGE
Formulation:	Lyophilized from a 0.2 μ m filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH7.4
Reconstitution:	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to a concentration of 0.1-1mg/ml by gently pipetting 2-3 times, don't vortex.
Storage:	The lyophilized protein is stable at -20°C for up to 1 year, the protein solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle.

Synonyms: Ankyrin repeat and SOCS box containing 13, Ankyrin repeat and SOCS box containing protein 13, Ankyrin repeat and SOCS box protein 13, Ankyrin repeat domain containing SOCS box protein Asb 13, Ankyrin repeat domain containing SOCS box protein Asb13, ASB 13, ASB-13, Asb13, ASB13_HUMAN, FLJ13134, MGC19879

SDS-PAGE:



Safety Note: This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.