
FineTest Lysis Buffer (for ELISA)

Product Information

Catalog No.:	E050
Size:	100mL
Form:	liquid
Storage:	Store at -20°C, valid for one year. When used frequently, it can be stored at 2-8°C for a short time.

Background

The FineTest Lysis Buffer produced is a traditional rapid lyse solution for tissue cells. Protein samples obtained by Finetest Lysis Buffer can be used for ELISA.

This product can be used on animal, plant cell or tissue samples. The main components are 50mM Tris(pH 7.4), 150mM NaCl, 0.1% SDS. When using this lysate, the user can add specific inhibitors or not according to the specific purpose.

Protein samples obtained by FineTest Lysis Buffer can use FineTest BCA protein assay kit (K001) to measure protein concentration. Due to the high concentration of detergent, the protein concentration of the sample obtained from the FineTest Lysis Buffer cannot be determined by Bradford method.

Note:

You need to bring your own PMSF. Add PMSF to a final concentration of 1mM a few minutes before use, or add the appropriate mixture of protease phosphatase inhibitors.

All steps for cracking samples are performed on ice or at 4°C.

This product is only used for scientific research by professionals, and shall not be used for clinical diagnosis or treatment, food or medicine, or stored in ordinary residences.

For your safety and health, please wear a lab coat and disposable gloves to operate.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com