
Poly HRP-conjugated Affinipure Goat Anti-Mouse IgG (H+L)

Product Information

Catalog No.:	ES-0003
Size:	100 μ l/500 μ l/1mL
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	HRP purified from horseradish root
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months.

Application

Tested Application:	ELISA, WB, IHC
Suggested dilution:	1:5000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates 1:50-1:400 for IHC

Note: This product is for research use only, not for diagnostic or therapeutic use.