

## Poly L-Lysine coated slides

### Product Information

Catalog No.: IHC-0016

Storage: Room temperature for 24 months.

Size: 50 slides /kit

Product Dimensions:  $76 \times 26 \times 1$ mm

### Introduction

The surface of this product is treated with Poly L-Lysine, the surface is positively charged, and it is easier to adhere to a variety of conventional slices or other appropriate samples.

### Application

Tested Application: IHC

Suggested dilution: Ready to use.

**Note: This product is for research use only, not for diagnostic or therapeutic use.**